

TECHNICAL DATA-SHEET

MOBILE PROTECTION

1 DOOR



P-PM-01

- Wall cm 80 x 190h
- Lead glass eq pb mm 2.0 dim. cm 24 x 30
- WEIGHT: 73 Kg

P-PM-02

- Wall cm 80 x 190h
- Lead glass eq pb mm 2.0 dim. cm 30 x 40
- WEIGHT: 75 Kg

3 DOOR



P-PM-03

- Wall cm 30 + 80 +30 x 190 h
- Lead glass eq pb mm 2.0 dim. cm 24 x 30
- WEIGHT: 105

P-PM-04

- Wall cm 30 + 80 +30 x 190 h
- Lead glass eq pb mm 2.0 dim. cm 30 X 40
- WEIGHT: 107 Kg

- Mobile protection, sliding on wheels equipped with brakes;
- Frame in tubular steel painted in focus;
- Antix panel Pb mm 2, coating with plastic laminate;
- Complete of lead glass eq pb mm 2.0.