

## TECHNICAL DATA-SHEET



- Homogeneous lighting by n°3 fluorescent tubes with colour temperature 6000° K.
- Pressing rolls film-holding system
- Stove enamelled steel structure
- Outside polished steel frame
- White monochromatic diffuser
- Equipped with terminal board for electrical connection.
- Predisposed for explosion-proof installation.
- Which this declaration refers, conforms to Directive 93/42/CEE – Medical Devices - Class 1  
Moreover meets the Directive 73/23 CEE (Low tension) 89/336 CEE (Electromagnetic Compatibility)

The construction conforms to Norms EN 60601 -1,2 - Din 6856 -1,2

## TECHNICAL DATA - SHEET

<b>CODE</b>	<b>SCHERMATURA</b>	<b>OUTSIDE DIMENSION</b>	<b>WEIGHT</b>
<b>S-SM-13</b>	Pb mm 1	CM 100 X 150 H	KG 56
<b>S-SM-14</b>	Pb mm 1	CM 120 X 150 H	KG 67
<b>S-SM-15</b>	Pb mm 2	CM 100 X 150 H	KG 74
<b>S-SM-16</b>	Pb mm 2	CM 120 X 150 H	KG 86