


Kalitenin tanımı... Definiton of quality



atrax

Protect Your Operating Theatre from Infections

By simply using different undercarriages inside and outside of the operating theatres, ATRAX TT helps you prevent the sterile area from being contaminated during the patient's transfer

Patient Security First

ATRAX Patient Transfer Stretcher TT is designed to provide for optimal safety and a smooth process during the transfer of the patient between sterile and nonsterile areas. The winning patented locking mechanism of ATRAX TT, you do not have to worry about the risk of patient falling or other discomforts. Our innovative lock ensures that both undercarriages are securely interlocked before letting the patient seamlessly slide into the next area.

General Specifications

- ** Side rails in circular down-swing motion
- ** Rising head section
- ** Secure brake system with patented locking mechanism
- ** One IV poles
- ** Flame retardant, viscoelastic water proof matters
- ** Side rails in circular swing-down motion
- ** Bumpers
- ** Transfusion pole



- ** Weight Capacity 240 kg
- ** Total Length : 4300 mm
- ** Length : 2150 mm
- ** Width : 780 mm
- ** Width of lying surface : 620 mm
- ** Length of Lying Surface : 1970 mm
- ** Mattress Area Width : 615 mm
- ** Mattress Area Depth : 100 mm
- ** Height (Without Mattress) : 83,5 mm
- ** Height (With Mattress) : 93,5 mm



Fabrika

1373. Sk. İvogsan Ağaç İşleri Sitesi
No: 30/1-2-3 Kat Yenimahalle / Ankara
Tlf: (312) 395 43 30 - Fks: (312) 395 17 59
www.atraxmedikal.com - info@atraxmedikal.com

Satış

Altay Mah. 2665 . Cad. No:1 E. Blok K:2 D: 17
Etimesgut / Ankara
Tlf: (312) 502 09 48 - Fks: (312) 502 09 49
www.atraxmedikal.com - satis@atraxmedikal.com